

# A Comprehensive Regulatory and Technical Blueprint for a Dual-Pathway Neuromorphic XR Therapy Platform

## Foundational Governance: The ISO 13485 Quality Management System and Hard-Coded Safety Invariants

The successful development and market entry of a high-impact neuromorphic XR therapy platform like DreamScape necessitates a governance framework that is both rigorous in its adherence to global medical device standards and innovative enough to accommodate novel technological pathways. The core of this framework is built upon two pillars: the universal adoption of the ISO 13485 standard for Quality Management Systems (QMS) and the proactive, multi-layered implementation of "hard-coded safety invariants," most notably the explicit exclusion of REM sleep integration. This dual approach ensures that while the platform explores cutting-edge applications in both clinical therapeutics and non-sleep research, it maintains an uncompromising commitment to patient safety, data integrity, and regulatory compliance throughout its entire lifecycle. The

DreamScape\_ResearchOrchestrator\_v1\_0\_0 ALN module serves as the formal mechanism for coordinating these parallel efforts, explicitly identifying the mandatory requirements for an ISO 13485-style QMS, ISO 14971 risk management, and stringent data privacy controls as the foundational principles governing all work <sup>46</sup> <sup>49</sup>. This strategic alignment is not merely a procedural exercise but a fundamental design philosophy that shapes every aspect of the system, from database schema to clinical protocols.

The bedrock of this governance structure is the ISO 13485:2016 standard, which has become the internationally recognized benchmark for QMS in the medical device industry <sup>1</sup> <sup>6</sup>. It provides a comprehensive framework covering all stages of the device lifecycle, from initial concept and design through production, servicing, and post-market surveillance <sup>1</sup>. Its importance is underscored by the U.S. Food and Drug Administration's (FDA) decision to harmonize its own Quality System Regulation (21 CFR Part 820) with ISO 13485, a move that will require U.S. companies to align their QMS with the international standard by February 2, 2026

2 42 . For a complex system like DreamScape, which involves software, hardware, and clinical data, adherence to ISO 13485 is non-negotiable. The standard mandates documented procedures for identification, storage, protection, retention, and disposition of records, ensuring that all critical information—from design files and manufacturing records to regulatory submissions—is meticulously controlled and traceable 2 43 . This directly supports the need for a centralized data lake and warehouse that de-identifies and stores therapy sessions, EEG data, and outcome measures, with calibration records tying each recording to specific devices and sites for cross-site comparability 3 7 . Furthermore, ISO 13485 requires the implementation of a risk management process conforming to ISO 14971:2019, making it the cornerstone of hazard identification and mitigation 13 51 . This creates a seamless integration between the overarching QMS and the specific risk management activities required for a Software as a Medical Device (SaMD) 49 .

A powerful and distinctive feature of the proposed architecture is the defense-in-depth strategy for safety, centered on the hard-coded prohibition against any form of REM sleep integration. This is not treated as a simple software filter but as a fundamental safety constraint embedded at every level of the system, reflecting a deep understanding of risk management principles outlined in ISO 14971, which prioritizes inherent safety-by-design over reliance on user warnings or information for safety alone 47 . At the database level, the SQL schema (`neuro_ar_lab_v1_0_0.sql`) enforces this invariant through a CHECK constraint on the `study` table, programmatically preventing any attempt to label a study or session as REM-related 17 47 . This is complemented by a corresponding behavior tree in the ALN policy module (`rem_exclusion_policy_v1_0_0.aln`), which codifies the same rule at the logical application layer, creating a redundant safeguard against accidental or malicious circumvention 47 . This architectural choice is critical because, under frameworks like the EU MDR and FDA guidance, SaMD that monitors physiological processes or impacts treatment decisions can be classified as Class IIa, IIb, or even III based on the severity of potential harm 19 20 . By proactively excluding REM sleep monitoring and intervention, the developers mitigate a vast array of potential hazards related to dream manipulation, seizure induction, or adverse psychological events associated with lucid dreaming or nightmare disorders, thereby lowering the overall risk profile of the platform 50 . This design choice is further reinforced at the hardware and protocol levels; the specifications for the AR gait lab hardware are explicitly configured for awake and ambulatory use only, and the clinical protocols focus exclusively on pre-sleep relaxation and sleep initiation dynamics rather than sleep architecture during REM

Calibration and traceability are other non-negotiable requirements deeply embedded within this governance framework. ISO 13485 Clause 7.6 explicitly mandates the control of monitoring and measuring equipment, requiring calibration at specified intervals or prior to use, with results traceable to national or international standards <sup>3</sup>. For DreamScape, this applies directly to the clinical-grade XR/BCI hardware, including EEG headsets, biosensors, and motion-capture systems <sup>6</sup>. The proposed database schema provides a robust blueprint for meeting this requirement, featuring dedicated tables for `device_asset`, `calibration_profile`, and `calibration_record`[[2,3]]. Each record must contain essential metadata such as equipment ID, calibration date, results, and the name of the person who performed the calibration, ensuring complete traceability for every measurement taken <sup>3</sup>. This meticulous record-keeping is vital for both clinical trials, where consistent data quality across multiple sites is paramount, and for post-market surveillance, where any anomalous data can be traced back to the specific device and its last calibration status <sup>70</sup>. The SOPs for the hospital gait lab, aligned with published protocols and ISO 13485 guidance, provide a practical example of how this is implemented, starting with a pre-check of device IDs and calibration due dates in the central database before any patient session can begin <sup>3</sup> <sup>4</sup>. This ensures that no device without a valid, in-date calibration record is used, forming a critical line of defense against inaccurate data collection and subsequent flawed clinical conclusions. The combination of a robust QMS, hard-coded safety invariants, and rigorous calibration procedures establishes a foundation of trust and reliability that is essential for navigating the complex landscape of digital health regulation.

Governance Component	ISO 13485 Requirement	Implementation in DreamScape Architecture
Quality Management System (QMS)	Establish, document, implement, and maintain a QMS. Ensure process-based continual improvement.	Adherence to ISO 13485:2016 is the guiding principle for all layers, from design controls to post-market surveillance. <a href="#">1</a> <a href="#">6</a> <a href="#">49</a>
Risk Management	Implement a process conforming to ISO 14971:2019 for identifying hazards, estimating risks, implementing controls, and evaluating residual risk.	Integrated throughout the development lifecycle, with a focus on inherent safety-by-design (e.g., REM exclusion). <a href="#">13</a> <a href="#">47</a> <a href="#">51</a>
Control of Monitoring Equipment	Calibrate monitoring equipment at specified intervals; ensure traceability of calibration status.	Centralized calibration_record table links EEG recordings to device IDs and calibration dates for cross-site comparability. <a href="#">2</a> <a href="#">3</a> <a href="#">6</a>
Document Control	Maintain documented procedures for the identification, storage, protection, retrieval, retention, and disposition of records.	Electronic Document Management System (EDMS) to manage all design files, specifications, and regulatory documents with version control. <a href="#">2</a> <a href="#">43</a>
Traceability	Ensure traceability of product throughout its lifecycle.	Calibration tables and audit logs tie each recording and event to device IDs, subject pseudonyms, and operator roles. <a href="#">2</a> <a href="#">3</a> <a href="#">41</a>
Hard-Coded Safety Invariant	Not explicitly stated, but implied by the need for inherent safety-by-design.	Multi-layered enforcement of REM sleep exclusion via SQL CHECK constraints, ALN policy logic, and hardware/software protocols. <a href="#">17</a> <a href="#">46</a> <a href="#">47</a>

## Phased Regulatory Pathways and Clinical Validation Strategies for Insomnia Treatment

The development of DreamScape is structured as a progressive journey along a clear regulatory pathway, transforming a wellness-focused virtual reality (VR) relaxation tool into a cleared prescription digital therapeutic (PDT) for chronic insomnia. This phased approach allows the company to build a robust technology stack, gather preliminary data, and lay the groundwork for increasingly stringent regulatory submissions. The roadmap begins with Phase 1, focused on evidence-aligned VR relaxation modules marketed as a wellness product, which likely falls under Class I FDA regulations [9](#) [10](#). While simpler from a regulatory perspective, this phase is strategically crucial for establishing a telemetry infrastructure, collecting usage metrics, and building a foundation of documentation that will be invaluable for future medical device positioning [21](#). Phase 2 marks the transition to a Software as a Medical Device (SaMD), introducing CBT-I inspired therapy sequences and a clinician dashboard [50](#). This step moves the product closer to a Class II classification, as it now performs a medical purpose—treating insomnia—

without being part of a physical hardware device, a definition adopted by regulators worldwide [8](#) [15](#) . The feasibility and pilot studies conducted in this phase prepare the ground for seeking regulatory clearance as an adjunctive insomnia treatment, leveraging the precedent set by products like Somryst® [55](#) [68](#) .

The true challenge and opportunity lie in Phases 3 and 4, where the introduction of EEG and closed-loop control represents a significant leap in risk and complexity. This evolution pushes the product towards a higher risk classification, potentially Class II or even Class III, depending on the specifics of the intervention and the potential consequences of failure [9](#) [18](#) . The ultimate goal of achieving regulatory-grade VR insomnia treatment status firmly places DreamScape in the category of a Prescription Digital Therapeutic (PDT), a sub-class of SaMD that delivers a therapeutic intervention and requires a physician's prescription for access [21](#) [63](#) . The existence of FDA-cleared DTx products such as Somryst® (formerly SHUTi) and SleepioRx provides a detailed and actionable roadmap for this process [58](#) [62](#) . These products were successfully cleared primarily through the 510(k) pathway, demonstrating substantial equivalence to a legally marketed predicate device and supported by extensive clinical trial data from randomized controlled trials (RCTs) involving thousands of patients [55](#) [67](#) [68](#) . This precedent validates the chosen path and provides a clear target for the type and volume of clinical evidence required for submission.

Central to this entire process is the three-pillar clinical evaluation framework established by the International Medical Device Regulators Forum (IMDRF) and adopted by the FDA [65](#) [72](#) . Any claim made by DreamScape, particularly those related to its neuromorphic EEG closed-loop features, must withstand scrutiny across three distinct domains. First, Valid Clinical Association requires establishing a scientific link between the technology's output (e.g., real-time arousal detection derived from EEG signals) and the targeted clinical condition of insomnia [72](#) . This involves a thorough literature review and alignment with established clinical constructs, such as sleep onset latency (SOL), wake after sleep onset (WASO), and total sleep time (TST) [54](#) . Second, Analytical Validation is the technical proof that the software correctly and reliably processes input data (the raw EEG signal) to generate accurate and precise outputs (the arousal probability score) [16](#) [64](#) . This phase would involve validating the performance of the neuromorphic classifier against a gold standard like polysomnography (PSG). The cited accuracy of 79% for the Sleeptracker-AI Monitor in epoch-by-epoch sleep staging provides a relevant benchmark for what can be expected from a non-clinical grade algorithm, highlighting the need for continuous refinement to meet the higher standards

required for a medical device 35 . Finally,Clinical Validationis the demonstration that the use of the technology leads to clinically meaningful outcomes in the target population 40 72 . This is the domain of traditional RCTs. The proposed validation study protocol for DreamScape is well-aligned with this framework, using the Insomnia Severity Index (ISI) as the primary endpoint, a measure shown to be effective in trials for existing DTx products 58 60 . Success in this final pillar is measured by quantifiable improvements in validated scales and adherence to benchmarks defined by clinical guidelines, such as the American Academy of Sleep Medicine (AASM), which recommends a clinically significant threshold of at least a 10-minute reduction in SOL 54 .

Phase	Primary Focus	Intended Use / Claim	Likely Regulatory Classification	Key Milestone
Phase 1	VR Relaxation	Wellness / Stress Reduction	Class I SaMD / Wellness Product 9 10	Launch of basic relaxation modules with self-reported outcome tracking.
Phase 2	Structured CBT-I Adjunct	Adjunct to Standard Care for Insomnia	Class II SaMD 9 50	Release of clinician dashboard and completion of pilot clinical studies.
Phase 3	EEG-Informed Closed Loop	Research Mode for Arousal Modulation	Class II SaMD (Research Use Only) 9	Collection of large datasets from research sites to refine algorithms.
Phase 4	Regulatory-Grade Treatment	Treatment of Chronic Insomnia	Class II Prescription Digital Therapeutic (PDT) 21 63	Submission for FDA 510(k) or De Novo clearance as a medical device.

The clinical validation strategy must also account for the nuances of insomnia treatment. EMA and FDA guidelines recommend trial durations of 4–8 weeks and often favor three-arm designs with placebo and active comparator arms 52 . The choice of control arm is critical; for insomnia-focused DTx, a "sham" control is often achieved by removing disease-specific CBT content while retaining engagement structures, a method accepted by the FDA 61 . Studies have shown that context-appropriate, non-therapeutic digital controls are feasible and accepted by regulators 58 . The primary outcome, typically the Insomnia Severity Index (ISI), is supplemented by secondary endpoints like sleep diary metrics (SOL, WASO), actigraphy-derived sleep measures, and validated anxiety/depression scales (e.g., GAD-7, PHQ-9) 44 57 . Adherence monitoring is also a key component, as seen in the Somryst protocol, which requires patients to complete a minimum number of sleep diaries to progress through treatment cores 57 60 . A robust statistical analysis plan, including intention-to-treat analysis and subgroup analyses, is essential to demonstrate efficacy conclusively 53 60 . By following this structured, evidence-based approach, DreamScape can navigate the complexities of SaMD regulation and establish itself as a safe and effective treatment for chronic insomnia.



# Deep Dive into Technical Architecture: From Neuromorphic Inference to Secure Data Platforms

The technical architecture of DreamScape is a sophisticated, multi-layered system designed not only to deliver therapeutic value but also to meet the stringent demands of a medical device environment. The hierarchical structure, comprising clinical knowledgebases, therapy orchestration, neuromorphic inference, and a robust data platform, enables a phased development and regulatory pathway [9](#) [21](#) . The most forward-looking and technically demanding layer is the neuromorphic inference engine, which is responsible for processing real-time EEG data to detect arousal states and enable a closed-loop adaptation of the VR experience [32](#) . This choice of technology is driven by the need for near-real-time, low-power inference directly on the headset or bedside unit, minimizing latency and power draw—a critical requirement for a wearable, always-on sensor system [33](#) [34](#) . Event-driven neuromorphic hardware or accelerators, such as the Intel Loihi 2 chip capable of simulating billions of neurons, are ideal substrates for this task [33](#) . Spiking Neural Networks (SNNs), the computational model underlying this hardware, process information through discrete, time-dependent spikes, allowing for asynchronous computation where periods of silence incur zero energy cost [34](#) . This native synergy between SNNs and neuromorphic chips makes them exceptionally well-suited for the sparse, event-driven nature of EEG signals, enabling continuous sleep-onset detection with ultra-low power consumption [34](#) .

However, the adoption of adaptive AI/ML presents unique challenges for regulatory compliance. To satisfy the traceability and documentation requirements of ISO 13485 and IEC 62304, the system must incorporate Good Machine Learning Practices (GMLP) and, critically, use explainable AI (XAI) models [13](#) [42](#) . An opaque "black box" model would fail to meet regulatory expectations, as it would be impossible to verify the reasoning behind its outputs or track changes effectively. Therefore, the development process must integrate XAI techniques to make the neuromorphic classifier's decision-making process transparent and auditable [42](#) . The closed-loop controller itself must operate within strict clinician-defined bounds, with a fallback mechanism to revert to conservative, simple relaxation content if the classifier confidence is low or unusual patterns are detected [32](#) . All closed-loop decisions must be logged for later clinical review and algorithm audits, providing a complete and immutable record of the system's behavior [17](#) . This logging capability is a core tenet of FDA guidance on electronic records, which emphasizes authenticity, integrity, and traceability throughout the full retention period [41](#) .

The data platform and analytics layer serves as the backbone of the entire system, supporting both product telemetry and clinical analysis <sup>30</sup>. It must be architected to handle a diverse range of data types, including de-identified therapy session logs, summaries of EEG and biosignals, sleep-stage trajectories, and validated outcome measures from sleep diaries and questionnaires <sup>7</sup>. A clinical data lake and warehouse is the appropriate architecture for this, providing a scalable and organized store for this rich dataset <sup>7</sup>. However, the most critical aspect of this layer is its implementation of privacy and security controls to comply with frameworks like HIPAA and GDPR <sup>37 71</sup>. Pseudonymization is a cornerstone of this strategy, defined as the replacement of direct identifiers while maintaining re-identification capability under controlled conditions <sup>7</sup>. This privacy-by-design measure is mandated by major privacy laws and must be implemented rigorously. Tools like WiseSpace or Mainzliste could be integrated to automate the de-identification process, applying configurable rules to remove or transform HIPAA identifiers and quasi-identifiers while classifying the re-identification risk <sup>7 39</sup>. This ensures that sensitive health information is protected, yet the data remains useful for research and analytics when accessed by authorized personnel under controlled conditions.

Furthermore, the data platform must facilitate continuous learning and change control, particularly for the AI/ML components. The FDA's guidance on Predetermined Change Control Plans (PCCPs) offers a regulatory pathway for updating AI/ML-enabled SaMD post-market without requiring a new submission for every modification <sup>30 71</sup>. Beacon Biosignals' SleepStageML, an FDA-cleared ML software for sleep staging, includes an authorized PCCP that covers four types of modifications, from retraining the model with new data to updating signal processing steps <sup>30 69</sup>. DreamScape's architecture should be designed to support a similar PCCP. This requires a robust CI/CD pipeline that can automatically trigger predefined verification and validation testing whenever a new algorithm version is deployed. The system must include version tracking for both the VR app and the classifiers, along with rollback capabilities to quickly revert to a stable version in case of issues <sup>70 71</sup>. This secure update mechanism, coupled with a locked "clinical" model for production use and experimental branches for research, allows for iterative improvement while maintaining regulatory integrity and ensuring patient safety <sup>71</sup>. The entire data infrastructure, from pseudonymization vaults to secure update servers, must be designed with these principles of privacy, traceability, and controlled change in mind to build a system that is not only clinically effective but also compliant and trustworthy.



Technical Component	Key Requirement	Implementation Strategy	Relevant Standards/Guidance
Neuromorphic Inference Engine	Low-latency, real-time, low-power EEG processing.	Use of event-driven neuromorphic hardware (e.g., Intel Loihi 2) and Spiking Neural Network (SNN) models.	N/A (Technological Choice)
Explainable AI (XAI)	Transparency and audibility of AI/ML model decisions.	Integrate XAI techniques to make classifier reasoning interpretable and traceable for regulatory review.	ISO 13485, IEC 62304, FDA GMLP <a href="#">13</a> <a href="#">42</a>
Closed-Loop Controller	Safe and effective adaptation of VR stimuli based on arousal.	Operate within clinician-defined bounds; implement a fallback to conservative content on low-confidence predictions.	ISO 14971 (Risk Management) <a href="#">32</a> <a href="#">47</a>
Data Storage	Scalable and organized storage of diverse clinical and telemetry data.	Clinical data lake/warehouse for storing de-identified session data, biosignals, and outcome measures.	HIPAA, GDPR <a href="#">7</a> <a href="#">37</a>
Privacy Protection	Compliance with data privacy laws for health information.	Mandatory implementation of robust pseudonymization to protect patient identity.	HIPAA, GDPR, EHDS <a href="#">7</a> <a href="#">39</a>
Change Control	Authorized and controlled updates to AI/ML algorithms post-market.	Development of a Predetermined Change Control Plan (PCCP) with automated testing and validation pipelines.	FDA PCCP Guidance, EU AI Act <a href="#">30</a> <a href="#">69</a> <a href="#">71</a>
Secure Updates	Provision for remote software updates with rollback capabilities.	Secure over-the-air (OTA) update mechanism with version control, encryption, and validation.	FDA Cybersecurity Guidance <a href="#">70</a> <a href="#">71</a>

## Hardware, Usability, and Post-Market Surveillance Integration

The successful translation of the DreamScape software and algorithmic stack into a safe and effective clinical product hinges on the careful integration of hardware, usability engineering, and post-market surveillance (PMS) planning from the earliest stages of development. The hardware specification for a clinical-grade AR/BCI headset must extend beyond raw performance metrics to encompass factors critical for long-term use, especially in a home environment. Key specifications include high-resolution displays with a sufficient field of view and low motion-to-photon latency to minimize cybersickness, which is a known hazard identified in ISO 14971 risk management processes [6](#) [50](#) . Inside-out tracking with six degrees of freedom (6-DoF) is essential for immersive experiences, and the inclusion of external cameras may be necessary for certain clinical tracking scenarios [6](#) . Ergonomics and hygiene are paramount; the headset must have adjustable fit and features like replaceable liners and cleanable surfaces to be suitable for multi-

patient use in clinics <sup>6</sup> . For the BCI component, an 8-32 channel dry or semi-dry EEG cap positioned to capture frontal and central activity is required for reliable sleep and arousal assessment <sup>6</sup> . This must be complemented by integrated sensors for photoplethysmography (PPG) and an inertial measurement unit (IMU) to refine sleep-stage detection through heart-rate variability (HRV) and movement metrics <sup>6</sup> . On-board low-power compute, potentially augmented by a neuromorphic accelerator, is necessary to run the inference models locally, reducing latency and dependency on constant connectivity .

Usability engineering is a mandatory component of the risk management process under standards like ISO 14971 and is a key element of the FDA's Quality System Regulation <sup>50 51</sup> . This process ensures that the device is not only effective but also safe and intuitive for the intended users—the patients—with minimal supervision from clinicians <sup>50</sup> . This involves conducting human-factors testing in realistic home environments to identify potential hazards related to misinterpretation of instructions, difficulty in setup, or adverse reactions like seizures or falls <sup>50</sup> . The design of the user interface, the clarity of the instructions provided to patients, and the effectiveness of the clinician dashboard for monitoring patient progress are all critical areas of focus. For example, the XR runtime must include safety overlays, such as a clear way to abort a session, and cybersickness monitoring tools that can prompt the user to take a break if necessary <sup>6</sup> . The entire system must be designed to be resilient, with robust physical design to minimize risks like cable entanglement, and electrical safety protocols aligned with medical device norms <sup>6</sup> . The feedback gathered during usability testing is fed directly back into the risk management file, allowing for iterative design improvements that enhance both safety and user experience <sup>51</sup> .

Post-market surveillance is not an afterthought but a core design requirement that must be planned for throughout the product lifecycle <sup>70</sup> . Effective PMS is mandated proactively by both EU MDR and FDA regulations and is essential for identifying unforeseen risks, monitoring long-term performance, and fulfilling obligations for continued safety and effectiveness <sup>51 70</sup> . The architecture of DreamScape must be built to facilitate PMS from day one. This includes integrating features for remote monitoring, which allows for the real-time collection of physiological and performance data, as well as automated error and event reporting <sup>70</sup> . Usage analytics are another critical component, providing insights into patient adherence, off-label use, and general system performance, all of which can be anonymized to protect patient privacy <sup>70</sup> . The secure update mechanism is perhaps the most

important PMS feature, as it allows for the rapid deployment of bug fixes, security patches, and minor improvements without compromising safety <sup>71</sup>. This must be managed through a robust change control process, with clear communication plans for users about what changes are being implemented <sup>69</sup>. The system should also be designed to autonomously collect real-world performance data, such as adverse events and end-user feedback, with minimal burden on the user, feeding this information back into the clinical evaluation and risk management processes <sup>40</sup>. This continuous learning loop ensures that the product remains safe and effective long after its initial launch, a key expectation for modern SaMD <sup>40</sup>. The ability to conduct real-world evidence generation is becoming increasingly important for regulatory authorities, and a well-designed PMS program is the foundation for generating this valuable data to support future regulatory submissions and maintain market authorization <sup>70</sup>.

## Strategic Synthesis: Coordinating Dual Architectures and Stakeholder Alignment

In conclusion, the provided materials outline a visionary and highly pragmatic blueprint for developing a next-generation neuromorphic XR therapy platform. The core strength of this strategy lies in its dual-pathway architecture, which simultaneously pursues two distinct markets under a unified governance structure. This approach allows the company to leverage its advanced technology in two ways: first, by developing theDreamScape sleep-therapy stack along a clear, defensible path toward FDA/EU MDR clearance as a prescription digital therapeutic for insomnia <sup>21</sup> <sup>63</sup>; and second, by establishing theREM-excluded neuromorphic AR gait lab stack as a high-impact research platform for the academic and rehabilitation community <sup>4</sup> <sup>6</sup>. This diversification is a brilliant strategic move, mitigating risk by not placing all development resources on a single regulatory path while creating multiple revenue streams. The gait lab can serve as a proving ground for the platform's neuromorphic and biosensing capabilities, generating valuable data and establishing credibility in a less-regulated space, while the DreamScape stack is being meticulously prepared for the rigors of medical device approval. The unifying thread that binds these two divergent paths is the ISO 13485-compliant QMS, which ensures that all development activities, regardless of the target use case, adhere to the highest standards of quality, safety, and traceability <sup>3</sup> <sup>49</sup>.

The success of this ambitious project depends on the disciplined execution of this dual-pathway strategy and a deep, early integration of regulatory thinking into the agile development process. Agile methodologies can be effectively employed for software development, but they must be adapted to accommodate the waterfall-style regulatory gates required for design controls and risk management <sup>46</sup>. The initial focus must be on establishing a robust QMS, creating a comprehensive Risk Management File per ISO 14971, and precisely defining the Intended Use, as this foundational document drives every subsequent decision regarding design, validation, and marketing claims <sup>47</sup>. The hard-coded safety invariants, particularly the multi-layered exclusion of REM sleep, should be viewed not as a limitation but as a core competitive advantage. This proactive, defense-in-depth approach to risk management demonstrates a profound commitment to patient safety that will resonate strongly with regulatory bodies and clinical stakeholders alike, potentially simplifying the classification and review process for the insomnia treatment application <sup>17 47</sup>.

Finally, effective communication tailored to the needs of different stakeholders is essential for aligning all teams toward a common goal. The DreamScape\_ResearchOrchestrator\_v1\_0\_0 module correctly identifies the necessity of this tailored approach <sup>46</sup>. For Regulatory Professionals, the conversation must center on the Quality Management System, the risk management file, and the clinical validation plan, demonstrating how the product meets all statutory and regulatory requirements. For Engineering Teams, the priority is the database schema, API definitions, CI/CD pipelines, and observability metrics that enable a compliant and scalable system. For Clinical Researchers, the emphasis must be on the study design, inclusion/exclusion criteria, validated outcome measures, and data quality protocols, assuring them that the platform is a reliable tool for generating meaningful clinical evidence <sup>46</sup>. For mixed stakeholder groups, a layered communication strategy—providing executive summaries, regulatory detail, technical annexes, and protocol synopses—is necessary to ensure everyone understands their role in achieving the shared objective. By maintaining an unwavering commitment to the foundational principles of ISO 13485, leveraging the hard-coded safety invariants as a strategic asset, and fostering clear, stakeholder-aligned communication, the project team can successfully navigate the complex journey from a promising concept to a world-class, clinically grounded neuromorphic XR therapy.

---

## Reference

1. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)
2. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)
3. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)
4. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)
5. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)
6. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWZya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdlYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZdZdaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)

wMDAtMDAwMC0wMDAwLXdIYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZEDaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA

7. [https://cdn.qwenlm.ai/qwen\\_url\\_parse\\_to\\_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWYya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdIYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZEDaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA](https://cdn.qwenlm.ai/qwen_url_parse_to_markdown/system00-0000-0000-0000-webUrlParser?key=eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJyZXNvdXJjZV91c2VyX2lkIjoicXdlb191cmxhcGFyc2VfdG9fbWYya2Rvd24iLCJyZXNvdXJjZV9pZCI6InN5c3RlbTAwLTAwMDAtMDAwMC0wMDAwLXdIYlVybFBhcnNlciIsInJlc291cmNlX2NoYXRfaWQiOm51bGx9.cz1eeZEDaQH5CgUaxwUmfEJfqTOZMoh3PbosHslSPA)
8. **Software as a Medical Device (SaMD)** <https://www.fda.gov/medical-devices/digital-health-center-excellence/software-medical-device-samd>
9. **Classify Your Medical Device** <https://www.fda.gov/medical-devices/overview-device-regulation/classify-your-medical-device>
10. **Software as a Medical Device FDA classification and ...** <https://qualysec.com/software-as-a-medical-device-fda/>
11. **Software as a Medical Device (SaMD) Documentation ...** <https://www.hardianhealth.com/insights/how-much-documentation-does-my-software-as-a-medical-device-really-need>
12. **SaMD Regulatory Pathways: A Step-by-Step Guide to ...** <https://sequenex.com/samd-regulatory-pathways/>
13. **Medical Device Risk Management for Software and SaMD | Enlil** <https://enlil.com/blog/medical-device-risk-management-for-software-and-samd/>
14. **SaMD: When Software is Considered as Device** <https://www.techmagic.co/blog/samd>
15. **An Introduction to Software as a Medical Device (SaMD)** <https://intuitionlabs.ai/articles/samd-definition-and-scope>
16. **Software as a Medical Device (SAMD) Clinical Evaluation** <https://topflightapps.com/ideas/samd-clinical-evaluation/>
17. **How to perform a clinical evaluation for Medical Device ...** <https://confinis.com/clinical-evaluation-for-medical-device-software/>
18. **SaMD: Software as a Medical Device** <https://orangesoft.co/blog/samd-software-as-a-medical-device>
19. **MDR Rule 11: Classification Rules for SaMD** <https://revolve.healthcare/blog/rule-11-mdr>
20. **Software As A Medical Device In Europe – New Regulatory ...** <https://www.jdsupra.com/legalnews/software-as-a-medical-device-in-europe-2712940/>
21. **Digital Therapeutics (DTx): Everything You Need to Know** <https://www.medidata.com/en/life-science-resources/medidata-blog/digital-therapeutics-dtx/>



22. Digital Therapeutics Market: 2025 Analysis & Growth Trends <https://intuitionlabs.ai/articles/digital-therapeutics-market-analysis>
23. Evaluating AI-driven characters in extended reality (XR) ... <https://www.sciencedirect.com/science/article/pii/S0933365725002052>
24. Physicians' Perceptions of AI and Extended Reality in ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC12607535/>
25. A data-centric artificial intelligent and extended reality ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC9434088/>
26. Software as a Medical Device (SaMD) <https://citrusbits.com/software-as-a-medical-device-samd/>
27. Is medical device regulatory compliance growing as fast ... <https://link.springer.com/article/10.1007/s12553-023-00775-x>
28. Evaluating the Efficacy and Safety of PROSOMNIA Sleep ... <https://clinicaltrials.gov/study/NCT06644573>
29. Approved trends and product characteristics of digital ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC12106666/>
30. Predetermined Change Control Plans: Guiding Principles for ... <https://ai.jmir.org/2025/1/e76854>
31. Machine Learning-Enabled Medical Devices Authorized by ... <https://www.mdpi.com/2227-9059/13/12/3005>
32. Generative Models of Brain Dynamics - PMC <https://pmc.ncbi.nlm.nih.gov/articles/PMC9335006/>
33. Spiking Neural Networks(SNN): Third-Generation AI and the ... <https://deepfa.ir/en/blog/spiking-neural-networks-neuromorphic-computing>
34. How Spiking Neural Networks Can Revolutionize Edge AI [https://www.linkedin.com/posts/shivajikant\\_tbeb-tbeb-spikingneuralnetworks-activity-7366450243007078400-4CII](https://www.linkedin.com/posts/shivajikant_tbeb-tbeb-spikingneuralnetworks-activity-7366450243007078400-4CII)
35. Polysomnographic validation of an under-mattress ... <https://www.sciencedirect.com/science/article/pii/S1389945722001368>
36. HIPAA and GDPR Compliance in IoT Healthcare Systems [https://www.researchgate.net/publication/379129933\\_HIPAA\\_and\\_GDPR\\_Compliance\\_in\\_IoT\\_Healthcare\\_Systems](https://www.researchgate.net/publication/379129933_HIPAA_and_GDPR_Compliance_in_IoT_Healthcare_Systems)
37. Privacy in electronic health records: a systematic mapping ... <https://link.springer.com/article/10.1007/s10389-022-01795-z>
38. Decentralized digital health ecosystems: a unified ... <https://www.frontiersin.org/journals/digital-health/articles/10.3389/fdgth.2025.1685628/full>

39. Possibilities for secondary data use of electronic health ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC12426038/>
40. SaMD Clinical Evaluation: What Do the FDA and IMDRF ... <https://www.greenlight.guru/blog/samd-clinical-evaluation>
41. FDA Guidelines for Managing Electronic Records <https://www.regdesk.co/blog/fda-qa-document-on-electronic-records-retention-and-management/>
42. ISO 13485 Implementation for Medical Device Companies <https://orangesoft.co/blog/iso-13485-standard>
43. Electronic Document Management Systems for Medical ... <https://www.cloudtheapp.com/electronic-document-management-systems-for-medical-devices/>
44. Digital cognitive-behavioural therapy application ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC11202642/>
45. Digital cognitive behavioral therapy for insomnia <https://aasm.org/digital-cognitive-behavioral-therapy-for-insomnia-platforms-and-characteristics/>
46. Using Existing Regulatory Frameworks to Apply Effective ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC10512990/>
47. The Risk Management + Design Controls Connection <https://www.greenlight.guru/blog/risk-management-design-controls>
48. Navigating regulatory challenges across the life cycle of a ... <https://www.sciencedirect.com/science/article/pii/S1532046425000851>
49. Design Controls: Requirements, Process, Flow Chart, ... <https://arterexmedical.com/medical-device-design-controls/>
50. Using Medical Device Standards for Design and Risk ... <https://biomedeng.jmir.org/2021/2/e26942>
51. Medical Device Risk Management: ISO 13485 & ISO 14971 <https://www.cognidox.com/blog/medical-device-risk-management-iso14971>
52. EMA and FDA psychiatric drug trial guidelines - PubMed Central <https://pmc.ncbi.nlm.nih.gov/articles/PMC8157504/>
53. Effect of Psychological and Medication Therapies for ... <https://jamanetwork.com/journals/jamanetworkopen/fullarticle/2813414>
54. Clinical Practice Guideline for the Pharmacologic Treatment of ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC5263087/>
55. Digital treatment for insomnia receives FDA clearance <https://aasm.org/digital-treatment-for-insomnia-receives-fda-clearance/>
56. Profile of Somryst Prescription Digital Therapeutic for ... <https://pubmed.ncbi.nlm.nih.gov/33226269/>

57. Somryst <https://dtxalliance.org/products/somryst/>
58. FDA-authorized software as a medical device in mental ... <https://www.nature.com/articles/s44184-025-00174-2>
59. Home | Somryst for Chronic Insomnia <https://www.somryst.com/>
60. Effect of a prescription digital therapeutic for chronic ... <https://www.frontiersin.org/journals/psychiatry/articles/10.3389/fpsy.2024.1450615/full>
61. Appropriate controls for digital therapeutic clinical trials <https://pmc.ncbi.nlm.nih.gov/articles/PMC9436387/>
62. The Effectiveness of Digital Cognitive Behavioral Therapy to ... <https://mental.jmir.org/2025/1/e84323>
63. Decoding FDA Labeling of Prescription Digital Therapeutics <https://pmc.ncbi.nlm.nih.gov/articles/PMC12090883/>
64. Key Insights from the FDA's DHT Guidance <https://www.medable.com/knowledge-center/blog-key-insights-from-the-fdas-dht-guidance>
65. FDA Issues Fourth and Final SaMD Clinical Evaluation ... <https://www.mwe.com/insights/fda-issues-fourth-final-software/>
66. Comparison of Product Features and Clinical Trial Designs ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC11530488/>
67. K191716 - 510(k) Premarket Notification - FDA <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm?ID=K191716>
68. In Pre-Cert 1st, FDA clears app-based treatment for chronic ... <https://www.medtechdive.com/news/pre-cert-program-fda-clears-pear-therapeutics-app-chronic-insomnia-digital-therapeutic/574965/>
69. Predetermined Change Control Plans: Guiding Principles ... <https://pmc.ncbi.nlm.nih.gov/articles/PMC12577744/>
70. How to Build Post-Market Surveillance into Device ... <https://sequenex.com/post-market-surveillance-software/>
71. AI Regulatory & Legal Frameworks for Biopharma in 2025 <https://intuitionlabs.ai/articles/ai-biopharma-regulation-2025>
72. 2017 FDA Guidance - Software as a Medical Device (SAMD) <https://innolitics.com/articles/fda-guidance-software-as-a-medical-device-samd-clinical-evaluation/>
73. Software as a Medical Device (SAMD): Clinical Evaluation <https://www.fda.gov/regulatory-information/search-fda-guidance-documents/software-medical-device-samd-clinical-evaluation>